

AIR FUMIGATION

Application with Chloritab as follows

Close Room fumigation in absence of Human presence

- Keep Chloritab tablet of 20g – take 500ml of water in plastic bowl and add 1 tablet .
- Keep it overnight .
- Small Aerator can be used for faster release .



Metal, Wood, Glass & Common Areas

- Spraying – 1g /10litres water – Machine used – Battery Operated sprayer/ hand sprayer/ Mist Blower .
- Through ULV FOGGER / Fogging Machine
- Keep Chloritab tablet of 5g – take 500ml of water in plastic bowl and add 1 tablet – Closed room Fumigation. (100SQ FT ROOM)

Fabrics, Sofas & Cloth

- Through ULV FOGGER / Fogging Machine
- Keep Chloritab tablet of 5g – take 500ml of water in plastic bowl and add 1 tablet – Closed room Fumigation. (100SQ FT ROOM)



Operation Theatres and Sensitive areas

- Through ULV FOGGER / Fogging Machine
- Keep Chloritab tablet of 20g – take 500ml of water in plastic bowl and add 1 tablet – Closed room Fumigation. (100SQ FT ROOM)

Street & open Areas

- Through ULV FOGGER / Fogging Machine
- Spraying – 10g /10litres water – Machine used - Battery Operated sprayer/ hand sprayer/ Mist Blower .



Cars, Buses ,Transport , Vehicles & Containers

- Through ULV FOGGER / Fogging Machine
- Spraying – 1g /10litres water – Machine used - Battery Operated sprayer/ hand sprayer/ Mist Blower

Global/ Domestic References & Recommendation

- Chinese Pharmaceutical Association – Expert Consensus on guidance and prevention strategies –Air Fumigation.
- US –EPA –United States Environment Protection Agency
- China National Health Commission.
- Gujarat State Dental Council
- National Institute of Virology Kerala Unit
- ICMR (Indian Council for Medical Research)